



PARK DENTAL SMILE OPTIMIZER INSTRUCTIONS FOR USE



Park Dental Smile Optimizer provides low frequency vibrations (30Hz) designed to help alleviate the pressure and discomfort associated with orthodontic braces and invisible aligners, while promoting faster results.

KIT CONTENTS:

- Smile Optimizer & Mouthpiece
- USB Cable
- Charger & Wall Socket Plug
- Charger Cradle

CHARGING INSTRUCTIONS:

- Place Smile Optimizer in charger cradle
- Insert USB cable into charger cradle
- Plug into power source
- Charge device for 2 hours
- The button will flash blue when a recharge is needed

OPERATION:

- Attach the included mouthpiece to the Smile Optimizer
- Press the device button to power on
- Place the device in your mouth using the recommended times below:
 - 10 minutes in the morning
 - 10 minutes in the evening

HOW TO CLEAN:

- Remove the mouthpiece from the device
- Clean mouthpiece with water and soft soap
- Rinse thoroughly

FAST RESULTS: Using the PD Smile Optimizer for just 20 minutes a day helps accelerate the movement of teeth and reduces the overall treatment time.

Visit us online for additional help, including tips and tutorials please visit

[Mysmile.ParkDentalAligners.com/Support](https://mysmile.parkdentalaligners.com/support)