



AORALSCAN 3

INTRAORAL SCANNER

+ LAPTOP
WITH SOFTWARE

\$13,999



- NVIDIA GTX2070 SUPER W/8GB RAM
- 512GB SSD | WINDOWS 10
- 15" ALIENWARE LAPTOP
- 6 CORE I7 PROCESSOR
- 32GB RAM



ENHANCE YOUR DIGITAL WORKFLOW

- ① AUTOMATIC ANTI-FOGGING
- ② ONE-BUTTON CONTROL
- ③ LIGHT WEIGHT (240g)
- ④ ERGONOMICALLY DESIGNED FOR OPTIMAL HAND-HELD OPERATION
- ⑤ DYNAMIC LED INDICATOR

25 SECONDS 30% FASTER

Single arch in 25 seconds.
Full mouth in 1 min.

AI SCAN SPEED UP THE SCAN

Automatically identifies and filters out unnecessary soft tissue data during scanning, resulting in a quicker and cleaner scan process.

2 TIP SIZES

Comfortably accommodates patients of all ages.

MOTION SENSING LESS CONTACT, MUCH SAFER

Dentists can scan with only limited contact with the computer, ensuring a safe and efficient scan experience.

30% MORE ACCURATE

Thanks to optimized algorithms and imaging mechanisms, Aoralscan 3 delivers results that are 30% more accurate, and better able to meet your clinical needs.



GET THIS SCANNER & LAPTOP FREE*

With the **Park Dental FlexScan Program**, you can add over **\$180,000 to your practice revenue while earning your own intraoral scanner and laptop at *no charge. Send us an average of 5-6 aligner cases per month and in approximately 18 months (100 cases total), the scanner and laptop are yours to keep and use as you please. Visit us online or contact us for more details.

*Scanner and laptop become the property of the clinician with purchase of 100 aligner cases.

**\$3,000 per case to the patient - \$1,195 (maximum) cost = \$1800 revenue per case



AORALSCAN 3

INTRAORAL SCANNER

SOFTWARE FEATURES:

- AUTOMATIC SCAN RETRACING
- UNIVERSAL APPLICATION
- MOTION SENSING
- AI SCAN



TECHNICAL SPECIFICATIONS

SCANNING FIELD	Standard scanner tip: 16 mm × 12 mm × 22 mm Mini scanner tip: 12 mm × 9 mm × 22 mm
SCAN PRINCIPLE	Non-contact scanner with structured light
DIMENSION (L x W x H)	281 mm × 33 mm × 46 mm
WEIGHT	240 ± 10 g (without cables)
OUT DATA	STL, OBJ, PLY
CONNECTOR	USB 3.0
POWER INPUT	12V DC/3 A

RECOMMENDED PC CONFIGURATION

CPU	Intel Core i7-8700 or higher
RAM	16 GB or more
HARD DRIVE	256 GB SSD or above
GRAPHIC CARD (GPU)	NVIDIA ® RTX 2060 6GB DDR3 or higher
OPERATING SYSTEM	Windows 10 Professional (64-bit) or later versions of Windows operating systems
SCREEN RESOLUTION	1920x1080, 60 Hz or higher
I/O PORTS	More than 2 Type-A USB 3.0 (or higher) ports



AoralScan is an intraoral scanner manufactured by Shining 3D and distributed by Park Dental Research Corp.

ORDER ONLINE AT
www.ParkDentalResearch.com/Shop or call:

1-800-243-7372

PM0227 Rev02